


Inverted Microscope

Venuslab CFM

High-end observation microscope dedicated to cell factories

- Ultra-long-distance contamination-free observation: Long-distance condenser, compatible with tall dishes/flasks, contamination-free
- HD optics & wide adaptability: Infinite flat-field objective + multi-spec stage, high-precision multi-scenario fit



One Platform Many Possibilities

Contact Us sales@venuslabtech.com

Get a Quote



Get Expert Advice
+65 8099 5547



Visit Us
www.venuslabtech.com

Overview

Cell factory observation: Professional microscopy solution

Venuslab CFM is an inverted microscope specifically designed for cell factory scenarios. Relying on a high-quality optical system and an adaptable structure, it not only meets the non-contamination observation needs of high culture containers but also is suitable for microscopic observation of cells/tissues in various fields such as scientific research and quality inspection, while balancing operational convenience and imaging reliability.

Features of Venuslab CFM

- 45° inclined binoculars for comfortable operation
- Adjustable interpupillary distance to reduce observation fatigue
- Supports the installation of phase contrast inserts for expandable functions
- Optional stage specifications to meet multiple needs
- High-rigidity anti-vibration structure, suitable for long-term observation

Service & Support

We are dedicated to delivering exceptional optoelectronic solutions to every client. From precision manufacturing and secure delivery to full-lifecycle technical support, we are here to ensure a seamless and reliable experience at every step.

1. Warranty Policy

Quality First, Worry-Free Operation

Warranty Period: We offer a **two-year** warranty service for all of our core optoelectronic products, effective from the date of shipment.

Coverage: We provide free repair or replacement services for malfunctions caused by material defects or workmanship errors under normal operating conditions.

Rapid Response: Upon receiving a warranty claim, we guarantee to initiate the assessment process within **24 hours** to minimize your equipment downtime.

2. Technical Support

Expert Team, Full-Process Guidance

Technical Consultation: Our team of senior optical engineers provides **24/7 online support** to assist with installation, commissioning, optical path alignment, and parameter optimization.

Scheduled Maintenance: We offer full-lifecycle maintenance recommendations, including firmware upgrades, optical component cleaning guidelines, and precision calibration services.

Training Services: We provide customized remote or on-site operational and safety training to ensure your team can operate the equipment efficiently and safely.

3. Logistics & Delivery

Precision Packaging, Global Reach

Professional Packaging: Given the fragile nature of optical instruments, we utilize industrial-grade shockproof, anti-static, and moisture-proof vacuum packaging to ensure zero damage during transit.

Logistics Partners: We partner with top-tier global logistics providers (**DHL / FedEx / UPS / SF Express**) to offer reliable shipping with real-time tracking.

Shipping Insurance: All shipments are fully insured to eliminate logistics risks.

4. Compliance & Certification


Strict Standards, Total Compliance

Quality Certification: Our manufacturing process is **ISO 9001 certified**, and our products comply with international standards such as **CE and RoHS**.

Export Compliance: "Committed to environmental responsibility, all our products comply with **RoHS 2.0 and REACH standards**, ensuring safety and global compliance."

Explore Series

Model	Glass circular stage plate size	Stage size
VL CFM-3	118mm	227mm X 208mm

 Get in touch with our team to explore configurations, request a quote, or learn more about customized solutions tailored to your needs.

Let us help you move science forward—faster and smarter.

[Get a Quote](#)



Get Expert Advice
+65 8099 5547



Visit Us
www.venuslabtech.com