

Electric Translation Stage

The VL Series High-Precision Electric Translation Stage is engineered for precise positioning applications in scientific research and industrial automation. Featuring a high-strength aluminum alloy body with a black anodized finish, combined with high-precision ball screws and linear guide systems, this series delivers exceptional linear motion accuracy, load capacity, and durability. Whether for optical alignment, semiconductor inspection, or biomedical instrumentation, the VL Series provides stable and reliable motion solutions.

Key Features

- **High-Precision Drive System**
Equipped with ground ball screws (or precision lead screws) and standard stepper motors for micron-level resolution and repeatability.
- **Superior Guide Mechanism**
Built-in precision linear guides (or cross-roller guides) minimize motion backlash, ensuring straightness and smooth operation under high loads.
- **Robust & Durable Construction**
Aviation-grade aluminum alloy with black anodized finish delivers corrosion and wear resistance for complex lab and industrial environments.
- **Flexible Expandability**
Standard mounting holes enable XY/XYZ multi-axis setups. Optional linear scales for full closed-loop control or servo motor upgrade for higher speeds.


Applications

- **Photonics & Laser Applications:** Beam alignment, fiber coupling, crystal adjustment.
- **Semiconductor & Microelectronics:** Wafer inspection, chip packaging positioning.
- **Biomedical:** Microscope stage movement, cell manipulation. **Automated Inspection:** Sensor calibration, custom automation equipment integration.

One Platform Many Possibilities

Contact Us sales@venuslabtech.com


Part 1: High Precision Motorized Linear Stage (VL-200UPTA Series)
Product Series: VL-200UPTA-(50-600)H



Model No.	Travel Range	Table Size	Load Capacity	Weight
VL-200UPTA-H50	50 mm	150×150 mm	20 kg	2.8 kg
VL-200UPTA-H100	100 mm	150×150 mm	20 kg	3.4 kg
VL-200UPTA-H150	150 mm	150×150 mm	30 kg	4.5 kg
VL-200UPTA-H200	200 mm	150×150 mm	30 kg	5.2 kg
VL-200UPTA-H300	300 mm	150×150 mm	30 kg	6.3 kg
VL-200UPTA-H400	400 mm	150×150 mm	50 kg	7.2 kg
VL-200UPTA-H500	500 mm	150×150 mm	50 kg	8.3 kg
VL-200UPTA-H600	600 mm	150×150 mm	50 kg	10.3 kg

Key Specifications:
Drive Mechanism: High Precision Ball Screw (Lead: 4mm/5mm)
Resolution: 1μm (20 Micro-steps) / 20μm (Full Step)
Repeatability: 1μm (Travel < 400mm) / 2μm (Travel ≥ 400mm)

Dimensions



Applications


- Precision Optical Alignment
- Laser Interferometry Systems
- Laboratory Sample Positioning

Part 2: Large Platform Motorized Linear Stage (VL-210UPTA Series)
Product Series: VL-210UPTA-(50-600)H



Model No.	Travel Range	Table Size	Load Capacity	Weight
VL-210UPTA-H50	50 mm	179×179 mm	20 kg	2.8 kg
VL-210UPTA-H100	100 mm	179×179 mm	20 kg	3.4 kg
VL-210UPTA-H150	150 mm	179×179 mm	30 kg	4.5 kg
VL-210UPTA-H200	200 mm	179×179 mm	30 kg	5.2 kg
VL-210UPTA-H300	300 mm	179×179 mm	30 kg	6.3 kg
VL-210UPTA-H400	400 mm	179×179 mm	50 kg	7.2 kg
VL-210UPTA-H500	500 mm	179×179 mm	50 kg	8.3 kg
VL-210UPTA-H600	600 mm	179×179 mm	50 kg	10.3 kg

Dimensions



Applications

- Semiconductor Wafer Inspection
- Bio-imaging & Scanning
- Multi-sensor Integration


Part 3: Compact Motorized Linear Stage (VL-220UPTA Series)

Product Series: VL-220UPTA-(50-300)H



Model No.	Travel Range	Table Size	Load Capacity	Weight
VL-220UPTA-H50	50 mm	115×115 mm	15 kg	2.8 kg
VL-220UPTA-H100	100 mm	115×115 mm	15 kg	3.4 kg
VL-220UPTA-H150	150 mm	115×115 mm	10 kg	4.5 kg
VL-220UPTA-H200	200 mm	115×115 mm	7 kg	5.2 kg
VL-220UPTA-H300	300 mm	115×115 mm	5 kg	6.3 kg


Dimensions



Applications

- Microscope XY Stages
- Fiber Coupling & Alignment
- Portable Device Integration

Part 4: Heavy Duty High Precision Linear Stage (VL-230HPTA Series)
Product Series: VL-230HPTA-(50-500)M




Model No.	Travel Range	Table Size	Load Capacity	Weight
VL-230HPTA-M50	50 mm	120×125 mm	20 kg	5μm
VL-230HPTA-M100	100 mm	120×125 mm	20 kg	5μm
VL-230HPTA-M150	150 mm	120×125 mm	30 kg	5μm
VL-230HPTA-M200	200 mm	120×125 mm	30 kg	5μm
VL-230HPTA-M300	300 mm	120×125 mm	30 kg	5μm
VL-230HPTA-M400	400 mm	160×160 mm	50 kg	5μm
VL-230HPTA-M500	500 mm	160×160 mm	50 kg	5μm
VL-230HPTA-M600	600 mm	160×160 mm	50 kg	5μm

Dimensions


Applications

VL-230HPTA-(50-300)M



- Industrial Automation Lines
- Large Object 3D Scanning
- Material Handling & Dispensing


VL-230HPTA-(400-500)M



Part 5: Specialized Series

1. Closed-Loop Series (Linear Scale Integrated)


Series Name: VL-GSxxTA (50-500) M Series



Description: Motorized stage equipped with an integrated linear scale for closed-loop control, achieving micron-level positioning accuracy.

Features: High position feedback accuracy, ideal for precision automation.

Dimensions



Applications


- Ultra-precision Metrology
- Lithography Process
- Closed-loop Feedback Control

2. Reverse Synchronous Dual Slider Series
Series Name: VL710TR150-SR



Drive Mode	Reverse Synchronous Dual Sliders
Travel Range	150 mm
Screw Type	Precision Ball Screw (Lead: 2mm)
Load Capacity	2 kg (Per Slider)

Dimensions



Applications

- Automatic Centering & Gripping
- Optical Slit Adjustment
- Symmetrical Material Testing

Part 6:Miniature Precision Motorized Stage (VL160TA Series)
Original Series: WN160TA
Feature: Compact 60x60mm platform, ideal for limited space applications.
Series Name: VL160TA (50-400)




Model No.	Travel Range	Table Size	Load Capacity	Weight
VL160TA50H	50 mm	60×60 mm	8 kg	1.4 kg
VL160TA100H	100 mm	60×60 mm	8 kg	1.7 kg
VL160TA150H	150 mm	60×60 mm	5 kg	2.3 kg
VL160TA200H	200 mm	60×60 mm	5 kg	2.6 kg
VL160TA300H	300 mm	60×60 mm	3 kg	3.0 kg
VL160TA400H	400 mm	60×60 mm	3 kg	4.5 kg


Key Specifications:

Drive Mechanism: Precision Ball Screw (Lead: 2mm)
Resolution: 10μm (Full step) / 0.5μm (1/20 Micro-step)
Max Speed: 10mm/sec
Guide Type: Precision Linear Guide


Dimensions




VL160TA50H




VL160TA100H



VL160TA150H



VL160TA250H



VL160TAH

Applications

- **Compact Optical Alignment:**
Ideal for adjusting lenses and mirrors in small optical setups.
- **Limited Space Positioning:**
Perfect for instruments where internal space is restricted.
- **Small Sample Scanning:**
Automated movement for small biological or material samples.

Part 7: Closed-Loop Stepper Motor Stage (VL180TA Series)
 Original Series: WN180TA
 Feature: 80x80mm platform equipped with a closed-loop stepper motor to prevent step loss.
Series Name: VL180TA (50-200)H




Model No.	Travel Range	Table Size	Load Capacity	Motor Type
VL180TA50H	50 mm	80×80 mm	15 kg	Closed-Loop Stepper
VL180TA100H	100 mm	80×80 mm	15 kg	Closed-Loop Stepper
VL180TA150H	150 mm	80×80 mm	10 kg	Closed-Loop Stepper
VL180TA200H	200 mm	80×80 mm	7 kg	Closed-Loop Stepper

Key Specifications:

Drive Mechanism: High Precision Ball Screw (Lead: 2mm)
 Repeatability: 1µm
 Max Speed: 12mm/sec


Dimensions



Applications


- **Reliable Automation:** Prevents step loss during long-duration automated tasks.
- **Feedback-Critical Positioning:** Used in processes where position verification is mandatory.
- **Vision Inspection Systems:** Ensures exact camera positioning without drift.

Part 8:High Load Cross-Roller Stage (VL260TH Series)
 Original Series: WN260TA
 Feature: Uses high-precision cross-roller guides for superior load capacity and stability.
Series Name: VL260TH Series



Model No.	Travel Range	Table Size	Load Capacity	Screw Lead
VL260TH30H	30 mm	80×80 mm	40 kg	1 mm
VL260TH50H	50 mm	120×125 mm	80 kg	4 mm
VL260TH100H	100 mm	120×165 mm	80 kg	4 mm
VL260TH150H	150 mm	165×245 mm	80 kg	4 mm

Dimensions



Applications

- **Heavy Load Material Handling:** Moving heavy industrial parts or fixtures (up to 80kg).
- **Industrial Production Lines:** Integration into assembly or dispensing machines.
- **Large Equipment Calibration:** Positioning heavy measurement devices.

Part 9:Product Series: VL261TH Series (Precision Compact)Original Series: VL261TA(15-75)H



Model No.	Travel Range	Table Size	Load Capacity	Screw Lead	Weight
VL261TH15H	15 mm	65 × 65 mm	10 kg	0.5 mm	1.2 kg
VL261TH30H	30 mm	90 × 90 mm	15 kg	1 mm	1.25 kg
VL261TH50H	50 mm	90 × 110 mm	15 kg	1 mm	1.4 kg
VL261TH75H	75 mm	90 × 110 mm	20 kg	1 mm	1.7 kg

Key Specifications:


Guide Type: Precision Cross-Roller Guide

Drive Mechanism: Precision Ground Screw

Repeatability: 3μm (for 15mm model) / 5μm (for others)

Resolution (20 Micro-step): 0.125μm (for 15mm model) / 0.25μm (for others)

Dimensions



Applications

- **Precision Optical Alignment :**
Ideal for fine-tuning lenses and mirrors requiring critical stability.
- **Micro-Assembly Operations :**
Stable movement for assembling small mechanical or electronic components.
- **Space-Constrained Inspection :**
High-accuracy positioning inside compact inspection instruments.

Part 10:Ultra-Compact Motorized Stage (VL-Micro TS)
 Original Series: WN262TA15H4040
 Feature: Extremely small footprint (40x40mm) for microscopic applications.
Series Name: VL-Micro TS




Model No.	Travel Range	Table Size	Load Capacity	Repeatability
VL-M262T15H4040	15 mm	40×40 mm	3 kg	1 μm

Key Specifications:

Drive: High Precision Ball Screw (Lead: 1mm)
 Guide: Cross Roller Guide

Dimensions



Applications


- **Microscopy Applications:** Precise sample movement under high-magnification microscopes.
- **Bio-chip Inspection:** Positioning alignment for microfluidic chips.
- **Fiber Optics:** Fine alignment of optical fibers in confined spaces.

Part 11: Configurable Precision Stage (VL-S262 Series)
 Original Series: WN262TA20...6060
 Feature: 60x60mm platform available in three precision grades (C/H/P).
Series Name: VL-S262-ECTA-20-6060



Model No.	Grade	Screw Type	Motor	Resolution
VL-S262-ECTA-20C-6060	Standard	Ball Screw	2-Phase Stepper	0.25μm
VL-S262-ECTA-20H-6060	Precision	Ground Ball Screw	2-Phase Stepper	0.25μm
VL-S262-ECTA-20P-6060	Ultra-Precision	Ground Ball Screw	5-Phase Stepper	0.1μm


Dimensions



Applications


- **General Laboratory Use (C-Grade):** Standard educational or experimental positioning.
- **Precision Optics (H-Grade):** Laser path construction and adjustment.
- **Semiconductor Inspection (P-Grade):** High-precision wafer or chip inspection requiring sub-micron accuracy.

Part 12:Rectangular Table Precision Stage (VL-262HPTA 6070 Series)Original Series: WN260TA30...6070
 Feature: 60x70mm rectangular table, available in C/H/P precision grades.
Series Name: VL-262HPTA-30-6070(C/H/P)



Model No.	Grade	Guide Class	Screw Class	Repeatability
VL-262HPTA-30C-6070	C-Grade	Cross Roller (C)	Ground Screw	5μm
VL-262HPTA-30H-6070	H-Grade	Cross Roller (H)	Ball Screw (C7)	3μm
VL-262HPTA-30P-6070	P-Grade	Cross Roller (P)	Ball Screw (C5)	1-2μm

Dimensions



Applications


- **Optical Slit Adjustment:** Adjusting light apertures or slits.
- **Sensor Calibration:** Moving sensors linearly for data collection.
- **Side-Mount Applications:** Suitable for mounting components that require a rectangular footprint.

Part 13:Low Profile Precision Stage (VL-262HPTA 8080 Series)
 Original Series: WN262TA...8080
 Feature: 80x80mm platform with a low-profile design using cross-roller guides.
Series Name: VL-262HPTA-(10-50)C-8080



Model No.	Travel Range	Table Size	Load Capacity	Weight
VL-262HPTA-10C-8080	10 mm	80×80 mm	10 kg	0.7 kg
VL-262HPTA-20C-8080	20 mm	80×80 mm	10 kg	0.7 kg
VL-262HPTA-30C-8080	30 mm	80×80 mm	10 kg	0.7 kg
VL-262HPTA-50C-8080	50 mm	80×80 mm	10 kg	0.7 kg


Dimensions



Applications


- **Multi-axis Stacking:** Ideal for building compact XY or XYZ systems due to its flat profile.
- **Height-Restricted Spaces:** Fits into machines where vertical clearance is very tight.
- **Stable Base Positioning:** Provides a wide, stable base for sensitive instruments.

Part 14:Center Drive Motorized Stage (VL-265CDTA Series)
 Original Series: WN265TA...
 Feature: Symmetrical structure with a central drive mechanism. Multiple sizes available.
Series Name: VL-265CDTA



Model No.	Travel	Table Size	Load Capacity	Screw Lead
VL-265CDTA-15H-4040	15 mm	40×40 mm	6 kgf	1 mm
VL-265CDTA-20H-6060	20 mm	60×60 mm	6 kgf	1 mm
VL-265CDTA-30H-8080	30 mm	80×80 mm	21 kgf	1 mm
VL-265CDTA-40H-100100	40 mm	100×100 mm	23 kgf	1 mm
VL-265CDTA-70H-120120	70 mm	120×120 mm	25 kgf	1 mm

Dimensions



Applications

- **Vertical (Z-axis) Applications:** Symmetrical design is excellent for vertical lifting stability.
- **Symmetrical Load Positioning:** Balanced movement for center-mounted loads.
- **Compact System Integration:** Central motor placement saves side space.

Part 15:Side Drive Thin Stage (VL-266SDTA Series)

Original Series: WN266TA...

Feature: Ultra-thin design achieved by placing the motor and drive mechanism on the side.


Series Name: VL-266SDTA




VL-266SDTA-20H-6060 VL-266SDTA-30H-8080 VL-266SDTA-40H-100100 VL-266SDTA-50H-120120

Model No.	Travel	Table Size	Load Capacity	Weight
VL-266SDTA-20H-6060	20 mm	60×60 mm	5 kgf	0.47 kg
VL-266SDTA-30H-8080	30 mm	80×80 mm	7 kgf	0.79 kg
VL-266SDTA-40H-100100	40 mm	100×100 mm	10 kgf	0.98 kg
VL-266SDTA-50H-120120	50 mm	120×120 mm	15 kgf	1.25 kg


Dimensions




VL-266SDTA-20H-6060



VL-266SDTA-30H-8080



VL-266SDTA-40H-100100




VL-266SDTA-50H-120120

Applications


- **Embedded System Integration:** The side-motor design allows the stage to be incredibly thin.
- **Microscope Sub-stages:** Fits easily between a microscope base and the sample holder.
- **Flat Panel Inspection:** Scanning flat surfaces where low stage height is critical.

Part 16:Nano-Adjustment Stage (VL-MiniTA Series)
Original Series: WN267TA20H
Feature: High-resolution stepper drive designed for nano-scale fine adjustments.
Series Name: VL-MiniTA Series



Model No.	Travel	Design Resolution	Repeatability	Load
VL-MiniTA-267-20H	20 mm	0.125 μm	< 1 μm	1 kg

Dimensions



Applications

- **Nano-scale Fine Tuning:** Ultra-high resolution for nanotechnology research.
- **Laser Interferometry:** Precise path length adjustment in interferometers.
- **Photonics Alignment:** Critical alignment of photonic integrated circuits.

 Get in touch with our team to explore configurations,request a quote,or learn more about customized solutions tailored to your needs.

Let us help you moves cience forward—faster ands marter.

Get a Quote

 Get Expert Advice
+65 8099 5547

 Visit Us
www.venuslabtech.com