

Objectives

M Plan Apo HR Series Objectives

High-order brightfield imaging + quantitative analysis

- Super Apochromatic Global Correction
- Synergy of High NA and Ultra-Long Working Distance
- Multi-Band High-Contrast Imaging
- Extreme Environment Adaptability Design



One Platform Many Possibilities

Contact Us sales@venuslabtech.com

Get a Quote



Get Expert Advice
+658099 5547 (WhatsApp)




Visit Us
www.venuslabtech.com

Overview

Introduction:

- Ultra-high resolution objective lens with an infinite conjugate parfocal distance of 95mm
- Chromatic aberration corrected within the visible spectrum range of 400-700nm
- Full field curvature correction for clear edge imaging
- Suitable for coaxial observation and convergence of visible lasers



Service & Support

We are dedicated to delivering exceptional optoelectronic solutions to every client. From precision manufacturing and secure delivery to full-lifecycle technical support, we are here to ensure a seamless and reliable experience at every step.

1. Warranty Policy

Quality First, Worry-Free Operation

Warranty Period: We offer a **two-year** warranty service for all of our core optoelectronic products, effective from the date of shipment.

Coverage: We provide free repair or replacement services for malfunctions caused by material defects or workmanship errors under normal operating conditions.

Rapid Response: Upon receiving a warranty claim, we guarantee to initiate the assessment process within **24 hours** to minimize your equipment downtime.

2. Technical Support

Expert Team, Full-Process Guidance

Technical Consultation: Our team of senior optical engineers provides **24/7 online support** to assist with installation, commissioning, optical path alignment, and parameter optimization.

Scheduled Maintenance: We offer full-lifecycle maintenance recommendations, including firmware upgrades, optical component cleaning guidelines, and precision calibration services.

Training Services: We provide customized remote or on-site operational and safety training to ensure your team can operate the equipment efficiently and safely.

3. Logistics & Delivery

Precision Packaging, Global Reach

Professional Packaging: Given the fragile nature of optical instruments, we utilize industrial-grade shockproof, anti-static, and moisture-proof vacuum packaging to ensure zero damage during transit.

Logistics Partners: We partner with top-tier global logistics providers (**DHL / FedEx / UPS / SF Express**) to offer reliable shipping with real-time tracking.

Shipping Insurance: All shipments are fully insured to eliminate logistics risks.

4. Compliance & Certification


Strict Standards, Total Compliance

Quality Certification: Our manufacturing process is **ISO 9001 certified**, and our products comply with international standards such as **CE** and **RoHS**.

Export Compliance: "Committed to environmental responsibility, all our products comply with **RoHS 2.0** and **REACH** standards, ensuring safety and global compliance."

Explore Series

Model	NA	\pm D.F	24 Eyepiece Field of View	WD
VL Plan Apo HR 5X	0.21	6.2 μ m	4.8mm	25mm
VL Plan Apo HR 10X	0.42	1.6 μ m	2.4mm	15mm
VL Plan Apo HR 20X	0.6	0.76 μ m	1.2mm	9.5mm
VL Plan Apo HR 50X	0.75	0.5 μ m	0.48mm	5.2mm
VL Plan Apo HR 100X	0.9	0.34 μ m	0.24mm	1.4mm

 Get in touch with our team to explore configurations, request a quote, or learn more about customized solutions tailored to your needs.

Let us help you move science forward—faster and smarter.

[Get a Quote](#)



Get Expert Advice
+658099 5547 (WhatsApp)



Visit Us
www.venuslabtech.com