

## Objectives

### Plan Fluor Ph Series Objectives

Widely used in live-cell research, cell analysis and sample fine imaging.

- Multimodal integration, reducing photobleaching
- High quality and efficiency in fluorescent imaging
- Outstanding details in phase contrast observation
- Wide sample compatibility



## One Platform Many Possibilities

Contact Us [sales@venuslabtech.com](mailto:sales@venuslabtech.com)

Get a Quote



Get Expert Advice  
+658099 5547 ( WhatsApp )



Visit Us  
[www.venuslabtech.com](http://www.venuslabtech.com)



### 3. Logistics & Delivery

#### Precision Packaging, Global Reach

**Professional Packaging:** Given the fragile nature of optical instruments, we utilize industrial-grade shockproof, anti-static, and moisture-proof vacuum packaging to ensure zero damage during transit.

**Logistics Partners:** We partner with top-tier global logistics providers (**DHL / FedEx / UPS / SF Express**) to offer reliable shipping with real-time tracking.

**Shipping Insurance:** All shipments are fully insured to eliminate logistics risks.

### 4. Compliance & Certification


#### Strict Standards, Total Compliance

**Quality Certification:** Our manufacturing process is **ISO 9001 certified**, and our products comply with international standards such as **CE** and **RoHS**.

**Export Compliance:** "Committed to environmental responsibility, all our products comply with **RoHS 2.0** and **REACH** standards, ensuring safety and global compliance."

### Explore Series

Model	NA	Focal Length	WD	Parfocal Distance
VL Plan Fluor PHL 4X	0.13	50mm	16.5mm	60.41
VL Plan Fluor PH1 10X	0.3	20mm	14.5mm	60.41
VL Plan Fluor PH1 20X	0.45	10mm	7.7-8.9mm	60.41
VL Plan Fluor PH2 40X	0.6	5mm	2.3-3.3mm	60.41
VL Plan Fluor PH2 60X	0.7	3.3mm	1.6-2.3mm	60.41

 Get in touch with our team to explore configurations, request a quote, or learn more about customized solutions tailored to your needs.

Let us help you move science forward—faster and smarter.

[Get a Quote](#)



Get Expert Advice  
+658099 5547 ( WhatsApp )



Visit Us  
[www.venuslabtech.com](http://www.venuslabtech.com)